

Elevate Sling

INSTRUCTIONS FOR USE

⚠ W A R N I N G S ⚠

- CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A LICENSED HEALTH CARE PRACTITIONER.
- CAUTION: FOR SINGLE PATIENT USE ONLY.
- CAUTION: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.
- WARNING: CAREFULLY READ FITTING INSTRUCTIONS AND WARNINGS PRIOR TO USE.
- WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.


Elevate LLC
(650) 862 2445
elevatebraces@gmail.com
025665 REV.0

elevate


www.elevatebraces.com

Fit the Elevate Brace


- 1 Place arm through shoulder loop and connect chest component (A) through D-Ring.



- 2 Adjust straps where shown to ensure the female magnetic attachment (B) is above the level of the heart.




- 3 The strap may be easily placed on either shoulder, depending on the comfort of the patient. The straps are cut to remove any excess material. The hook and loop fasteners are removed before cutting and reattached after the straps are cut to size.




Wrist Portion

- 4 Wrap the wrist band around the injured wrist with the male magnetic component (C) centered on the palm-side of the wrist.



Daily Use

Attach the wrist component to the shoulder strap by bringing the male magnetic component (C) straight-on to female attachment point (B). A distinct "click" confirms attachment. Pull straight up to release.



Washing Instructions: Fasten hook and loop closures. Hand wash in warm water using a mild detergent. Air dry.